




# Bandages and orthoses

for children and adults

Catalogue 2016 / 2017



# CONTENTS

**kids** Adults



## NECK

Cervicare .....	6	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Cervicare plus.....	7	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Cervicare stabil without trachea.....	8	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Cervicare stabil with trachea .....	9	<input type="checkbox"/>	<input checked="" type="checkbox"/>



## SHOULDER

Claviflex anatomic .....	12	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Claviflex straight.....	13	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Deltafix .....	14	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Gilchrist support .....	15	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Deltacare plus.....	16	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Deltacare women.....	17	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Deltacare Kissen.....	18	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Deltaflex.....	19	<input type="checkbox"/>	<input checked="" type="checkbox"/>



## ELBOW

Epicare plus.....	22	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Ellcare classic.....	23	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Ellcare flex .....	24	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Epicare .....	25	<input type="checkbox"/>	<input checked="" type="checkbox"/>



## HAND

Carpfix plus.....	28	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Manucare plus .....	29	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Carpus classic plus .....	30	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Carpus stabil plus .....	31	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Polli classic .....	32	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Polli stabil plus.....	33	<input type="checkbox"/>	<input checked="" type="checkbox"/>



## FINGER

Pollifix .....	36	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Fingerfix.....	37	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Stack finger splint .....	38	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Stack finger splint - Starter-Kit.....	38	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Stack splint for button hole deformity.....	38	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Stack splint for button hole deformity – Starter Kit	38	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Finger splint U-shape.....	39	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Frog finger splint.....	40	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Baseball finger splint .....	41	<input type="checkbox"/>	<input checked="" type="checkbox"/>



**SPINE**

Lumbocare basic.....	44	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Lumbocare basic plus .....	45	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>



**KNEE**

Inocare knee splint .....	48	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Helpcare .....	49	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Helpcare super .....	50	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Provicare.....	51	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Provicare plus.....	52	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Provicare senior .....	53	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Genocare plus .....	54	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Genocare ROM .....	55	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Fixcare 0°.....	56	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Fixcare 10°.....	56	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Fixcare 20°.....	56	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Patellacare .....	57	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Fourcare basic.....	58	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Fourcare classic.....	59	<input type="checkbox"/>	<input checked="" type="checkbox"/>



**LOWER LEG**

Kids' Walker .....	62	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Walker classic / Walker classic low .....	63	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Walker One Piece.....	64	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Walker ROM.....	65	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Walker air.....	66	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Walker ROM air.....	67	<input type="checkbox"/>	<input checked="" type="checkbox"/>



**ANKLE JOINT**

Castcare .....	70	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Malleocare plus .....	71	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>



**FOOT**

Universal kids' sandal .....	74	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Universal kids's shoe .....	74	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Forefoot relief shoe.....	75	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Universal-/Post surgery shoe .....	76	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Plaster cast shoe.....	77	<input type="checkbox"/>	<input checked="" type="checkbox"/>









# NECK

**kids** Adults

	<b>Cervicare</b> Anatomically formed cervical support with velcro fastener <b>p. 6</b>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
---	--	-------------------------------------	-------------------------------------

	<b>Cervicare plus</b> Anatomically formed cervical support with velcro fastener, reinforced design <b>p. 7</b>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
---	--	-------------------------------------	-------------------------------------

	<b>Cervicare stabil without trachea</b> Anatomically formed cervical support with enhanced stabilisation without trachea opening <b>p. 8</b>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
---	--	-------------------------------------	--------------------------

	<b>Cervicare stabil with trachea</b> Anatomically formed cervical support with enhanced stabilisation with trachea opening <b>p. 9</b>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
---	--	--------------------------	-------------------------------------

# Cervicare

Item. no. 111

## Anatomically formed cervical support with velcro fastener

The cervical support has a soft foam core and a fixed cover made of 100% cotton making it very skin friendly.

Very comfortable to wear due to cavity on the chin.  
Easily put on with velcro fastener.

Colours as shown on the picture. Children's sizes with velcro strap in different colours (non eligible).



Children



Adults

### Indication

- Traumatic disorders
- Acute or chronic cervical syndrome
- Trauma of cervical spine without bony or discoligamentary involvement
- Cervical spine disorder of degenerative origin
- Acute muscular torticollis

### Effectiveness

- Limits reclination and tilt to one side
- Cervical spine relief

### Measuring

**Look straight**

**U = neck circumference**

**H = height of chin**

(from chin to sternoclavicular joint)



### Materials

#### Foam with cotton covering

Core made of soft foam  
Cover made of 100% cotton  
skin friendly  
latex free

### Size table Cervicare



Adults

Size	XXXS	XXS	XS	S	M	L	XL	XXL	XXXL
Circum. (U) in cm	21 - 25	25 - 29	29 - 33	33 - 36	36 - 39	39 - 42	42 - 46	46 - 50	50 - 54
Height (H) in cm	4 or 6			8 or 10					

# Cervicare plus

Item. no. 112

## Anatomically formed cervical support with velcro fastener, reinforced design

The cervical support has a soft foam core and a fixed cover made of 100% cotton making it very skin friendly.

Reinforced by an additional seam . Very comfortable to wear due to cavity on the chin. Easily put on with velcro fastener.

Colours as shown on the picture. Children's sizes with velcro strap in different colours (non eligible).



Children



Adults

### Indication

- Traumatic disorders
- Acute or chronic cervical syndrome
- Trauma of cervical spine without bony or discoligamentary involvement
- Cervical spine disorder of degenerative origin
- Rheumatic torticollis
- Acute muscular torticollis
- After cervical spine surgery

### Effectiveness

- Restricted mobility
- Limits reclination and tilt to one side
- Cervical spine relief and stabilisation

### Measuring



**Look straight**

**U = neck circumference**  
**H = height of chin**

(from chin to sternoclavicular joint)

### Materials

#### Foam with cotton covering

Core made of soft foam  
 Cover made of 100% cotton  
 skin friendly  
 latex free

### Size table Cervicare plus



Adults

Size	XXXS	XXS	XS	S	M	L	XL	XXL	XXXL
<b>Circum. (U) in cm</b>	21 - 25	25 - 29	29 - 33	33 - 36	36 - 39	39 - 42	42 - 46	46 - 50	50 - 54
<b>Height (H) in cm</b>	4 or 6			8 or 10					



# Cervicare stabil without trachea

Item. no. 122

## Anatomically formed cervical support with enhanced stabilisation via two inflexible synthetic half shells

Chin support and occipital support for comfortable and secure stabilisation. This cervical support has no trachea opening.

Colour as shown in the picture.

Cotton covers available as accessories.



Children


### Indication

- Suspected bony cervical injury
- After cervical spine surgery

### Effectiveness

- Total fixation of the cervical spine

### Measuring



#### Look straight

**U = neck circumference**

**H = height of chin**

**chin**

(from chin to sternoclavicular joint)

### Materials

#### Supporting shell material

stabilizing and adaptable  
skin friendly  
easy to clean  
latex free

### Size table

Cervicare stabil w/o trachea



Size	XXS	XS
Circumf. (U) in cm	23 - 26	26 - 29
Height (H) in cm	5 or 6.5	

# Cervicare stabil with trachea

Item. no. 121

## Anatomically formed cervical support with enhanced stabilisation with trachea opening

Two inflexible synthetic half shells offer chin support and occipital support for comfortable and secure stabilisation.

Colour as shown in the picture.

Cotton covers available as accessories.



Adults


### Indication

- Suspected bony cervical injury
- After cervical spine surgery

### Effectiveness

- Total fixation of the cervical spine

### Measuring



**Look straight**  
**U = neck circumference**  
**H = height of chin**  
 (from chin to sternoclavicular joint)

### Materials

**Supporting shell material**  
 stabilizing and adaptable  
 skin friendly  
 easy to clean  
 latex free

### Size table

Cervicare stabil with trachea Adults






Size	S	M	L	XL
<b>Circum. (U) in cm</b>	25 - 31	31 - 38	38 - 45	45 - 53
<b>Height (H) in cm</b>	6, 8, 11 or 13			







Adults

	<p><b>Claviflex anatomic</b> Anatomic clavicle support p. 12</p>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
	<p><b>Claviflex straight</b> Straight clavicle support p. 13</p>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Deltafix</b> Highly immobilising shoulder vest p. 14</p>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Gilchrist support</b> Shoulder and arm orthosis p. 15</p>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Deltacare plus</b> Shoulder-arm-adduction support with forearm pouch openable towards the front p. 16</p>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Deltacare women</b> Shoulder-arm adduction support with rotatable upper arm sling p. 17</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Deltacare Kissen</b> Waist belt with anatomically formed waist cushion, 15° abduction wedge and thumb hole p. 18</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Deltaflex</b> Stable shoulder-arm abduction support for immobilization of the shoulder joint p. 19</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>

# Claviflex anatomic

Item. no. 232

## Anatomic clavicle support

Backpack like support for stabilization and immobilization of the shoulder position after clavicle fractures.

The straps in the anatomic variety are especially designed so they are more comfortable to wear by children.

The straps can be adjusted at four points.

Colours as shown in the picture. Children's sizes with velcro strap in different colours (non-eligible).



Children


### Indication

- Clavicle fractures

### Effectiveness

Fixation after fracture  
Posture correction

### Measuring



**U =**  
**circumference of**  
**the chest**

Please measure as  
shown in the  
picture

### Materials

**3 layered cotton blend fabric**  
high density  
very stable  
hard-wearing  
especially skin friendly  
warming  
latex free

### Size table

#### Claviflex anatomic



Size	XXXS	XXS	XS
Circumf. (U) in cm	< 40	40 - 50	50 - 70

# Claviflex straight

Item. no. 231

## Straight clavicle support

Backpack like support for stabilization and immobilization of the shoulder position after clavicle fractures.

The straps can be adjusted at two points.

Colours as shown in the picture. Children's sizes with velcro strap in different colours (non-eligible).



Children



Adults


### Indication

- Clavicle fractures

### Effectiveness

- Fixation after fracture
- Posture correction

### Measuring



**U = circumference of the chest**

Please measure as shown in the picture

### Materials

**3 layered cotton blend fabric**  
 high density very stable hard-wearing especially skin friendly warming latex free

### Size table

#### Claviflex straight



Adults

Size	XXXS	XXS	XS	S	M	L	XL	XXL
Circumf. (U) in cm	< 40	40 - 50	50 - 70	50 - 70	70 - 105	105 - 125	125 - 140	140 - 155



# Deltafix

Item no. 221

**Highly immobilising shoulder vest made of strong cotton material**

Children's sizes with velcro strap in different colours (non-eligible). Adult's sizes as shown in the picture.

How to put it on:



**kids** + Adults



## Indication

- Subcapital fracture of the humerus
- Immobilisation of the shoulder after surgery
- After luxation and luxation fractures
- Muscle and tendon injuries (upper arm)
- Immobilisation of extra articular irritations and illnesses in connection with the Gleno-Humeral joint
- Bursitis acromialis
- Tendinitis of the long biceps tendon
- After surgery for reconstruction of the rotator cuff
- After arthroscopic removal of rotator cuff calcifications
- Conditions after AC joint ruptures
- Supraspinatus syndrome
- Defilee extensions (acc. to Neer) in impingement syndrome


## Effectiveness

- High fixation or immobilisation of shoulder and elbow

## Measuring

**U = circumference of the chest**

Please measure as shown in the picture



## Materials

### 3 layered cotton blend fabric

high density very stable hard-wearing especially skin friendly warming latex free

## Size table Deltafix

**kids**

Adults

Size	XXS	XS	S	M	L
<b>Circum. (U) in cm</b>	45 - 55	55 - 70	70 - 85	85 - 105	105 - 125
<b>Side</b>	left / right		left / right		

# Gilchrist support

Item no. 241



## Shoulder and arm orthosis

For immobilisation of arm and shoulder at the same time.

The angle of the arm can be adjusted through flexible velcro fasteners. Hand can move freely, even with immobilised elbow.

Children's sizes with velcro strap in different colours (non-eligible). Adult's sizes as shown in the picture.


### Indication

- Injuries of the acromioclavicular joint
- Soft tissue irritation in the shoulder area
- Muscle and tendon injuries as well as fractures of the upper arm
- Fracture of the shoulder blade
- Tendinitis of the long biceps tendon
- Supraspinatus syndrome
- Bursitis acromialis
- Luxation or luxation fracture
- Subcapital fracture of the humerus
- Immobilisation of extra articular irritations and illnesses in connection with the Gleno-Humeral joint
- Calcification of rotator cuffs
- Regenerative support after surgery

### Effectiveness

- Fixation or immobilisation of shoulder and elbow

### Measuring



### Materials

**3 layered cotton blend fabric**  
 high density very stable hard-wearing especially skin friendly warming latex free

Size table  
 Gilchrist support



Adults

Size	XXS	XS	S	M	L	XL	XXL
Length (L) in cm	30 - 38	33 - 44	45 - 52	52 - 66	54 - 70	56 - 76	56 - 78
Circumf. (U) in cm	14 - 28	14 - 28	25 - 32	28 - 36	34 - 40	38 - 48	42 - 54

# Deltacare plus

Item no. 212

## Shoulder-arm-adduction support with forearm pouch openable towards the front

Shoulder-arm adduction support for immobilization of the shoulder area. Easy-to-open forearm pouch with velcro fastener. This feature enables the patient to put on the support without any further help. The angle of the arm can be adjusted through flexible velcro fasteners.

Colours as shown on the picture.  
Children's sizes with velcro strap in different colours (non eligible).



**kids** + Adults



### Indication

- Injuries of the acromioclavicular joint
- Soft tissue irritation in the shoulder area
- Muscle and tendon injuries of the upper arm
- Immobilisation of extra articular irritations and illnesses in connection with the Gleno-Humeral joint
- Supraspinatus syndrome
- Bursitis acromialis
- Tendinitis of the long biceps tendon
- Calcification of rotator cuffs


### Effectiveness

- Fixation or immobilisation of shoulder and elbow

### Measuring

**U = circumference of the chest**

Please measure as shown in the picture



### Materials

**3 layered cotton blend fabric**  
high density very stable hard-wearing especially skin friendly warming latex free

Size table  
Deltacare plus

**kids**

Adults

Size	XXS	XS	S	UNI	XL
Circumf. (U) in cm	< 50	50 - 70	60 - 85	85 - 120	> 120

# Deltacare women

Item no. 213

## Shoulder-arm adduction support with rotatable upper arm sling

Shoulder-arm adduction support for immobilization of the shoulder area. Easy-to-open forearm pouch with velcro fastener. This feature enables the patient to put on the support without any further help. The angle of the arm can be adjusted through flexible velcro fasteners.

The upper arm sling is rotatable and can be moved from the upper arm to the forearm. This can be especially helpful for women with a larger circumference of the chest.

Colour as shown in the picture.



Adults

### Indication

- Injuries of the acromioclavicular joint
- Soft tissue irritation in the shoulder area
- Muscle and tendon injuries of the upper arm
- Immobilisation of extra articular irritations and illnesses in connection with the Gleno-Humeral joint
- Supraspinatus syndrome
- Bursitis acromialis
- Tendinitis of the long biceps tendon
- Calcification of rotator cuffs


### Effectiveness

- Fixation or immobilisation of shoulder and elbow

### Measuring

**U = circumference of the chest**

Please measure as shown in the picture



### Materials

**3 layered cotton blend fabric**  
 high density very stable hard-wearing especially skin friendly warming latex free

### Size table

#### Deltacare women

Adults

Size	S	UNI	XL
Circumf. (U) in cm	60 - 85	85 - 120	> 120

# Deltacare Kissen

Item no. 218

**Waist belt with anatomically formed waist cushion, 15° abduction wedge and thumb hole**

Comfortable to wear due to padded shoulder strap. Flexible velcro fastener.

Colour as shown in the picture.



Adults


## Indication

- Immobilisation of the shoulder after surgery
- Muscle and tendon injuries (upper arm)
- Immobilisation of extra articular irritations and illnesses in connection with the Gleno-Humeral joint
- After surgery for reconstruction of the rotator cuff
- Postoperative immobilisation After arthroscopic removal of rotator cuff calcifications

## Effectiveness

- Relief of ventral limbus and capsule

## Measuring



**L = Length of forearm**  
Measure arm length as shown

## Materials

**3 layered cotton blend fabric**  
high density very  
stable hard-wearing  
especially skin friendly  
warming  
latex free

## Size table

### Deltacare Kissen

Adults

Size	S	M	L	XL
Length (L) in cm	25 - 32	32 - 39	39 - 46	46 - 50



# Deltaflex

Item no. 219

## Stable shoulder-arm abduction support for immobilization of the shoulder joint

This support can be adjusted to six different abduction positions: 0°, 30°, 45°, 60°, 75° and 90°. Can be used left and right by folding the arm support to the other side and attaching the forearm and wrist strips to the other side. A waist and a shoulder strap provide for the necessary stability.

Colour as shown in the picture.



Adults


### Indication

- Immobilisation of the shoulder after surgery
- Muscle and tendon injuries (upper arm)
- Immobilisation of extra articular irritations and illnesses in connection with the Gleno-Humeral joint
- After surgery for reconstruction of the rotator cuff
- Postoperative immobilisation after arthroscopic removal of rotator cuff calcifications
- Shoulder dislocation
- After shoulder implants
- Subcapital fracture of the humerus
- Rupture of the acromioclavicular joint
- SLAP lesion

### Effectiveness

- Relief of ventral limbus and capsule
- Postoperative or posttraumatic relief of the shoulder
- Relief of tendons and ligaments

### Measuring



**L = Length of arm**  
Measure from shoulder joint to MCP joints

### Materials

**3 layered cotton blend fabric**  
high density very stable hard-wearing especially skin friendly warming latex free

### Size table Deltaflex

Adults

Size	S	L
Length (L) in cm	48 - 58	58 - 68



# ELBOW



Adults



## Epicare plus

Active support for elbow compression with 2 incorporated pads

p. 22



## Ellicare classic

Elbow orthosis with pinjoint

p. 23



## Ellicare flex

Elbow orthosis with incorporated aluminium splint, adjustable

p. 24



## Epicare

Epicondylitis brace

p. 25



## Epicare plus

Item no. 312

### Active support for elbow compression with 2 incorporated pads

This active support helps to reduce pain and release tensions. Every move leads to a massaging effect through the pads. This activates the muscles and speeds up the healing. The elbow joint can still be moved freely.

Excellent wear with highly elastic fabric, no slipping

Colour as shown in the picture.




### Indication

- Epicondylitis (chronic or acute)
- Repetitive Strain Injury (RSI)
- Tennis elbow, golfers' elbow, mouse elbow
- Arthrosis / arthritis
- Irritation after surgery
- posttraumatic
- Fibromyalgia

### Effectiveness

- pain-relieving
- tension-relieving
- massaging
- supports healing

### Measuring



**U = Circumference elbow**  
Measure as shown

### Materials

**3 layered synthetic blend fabric**  
breathable  
elastic skin  
friendly easy  
to clean  
Pads: 100% silicon

Size table  
Epicare plus

kids

Adults

Size	XXS	XS	S	M	L	XL	XXL
Circumf. (U) in cm	17 - 19	19 - 21	21 - 23	23 - 25	25 - 27	27 - 29	29 - 32

# Ellcare classic

Item no. 321

## Elbow orthosis with pin joint

Elbow orthosis with small pin joint for adjusting extension and flexion in steps of 10 degrees:

- Extension: 0° - 90°
- Flexion: 0° - 120°

Easily adjustable through velcro straps, even in swollen conditions.

Colours as shown in the picture. Children's sizes with velcro strap in different colours (non-eligible).



Children



Adults

## Indication

- Bursitis
- Olecranon trauma
- After surgery
- Instability of the elbow joint
- Posttraumatic and postoperative soft tissue injuries in and around the elbow
- Joint reconstruction
- Joint effusion and swelling

## Effectiveness

- Immobilisation of the elbow joint
- high stabilisation of side movements

## Measuring



**U = Circumference upper arm**

Measure as shown

**L = overall length of orthosis**

In relation to the overall length the joint is in the middle.

## Materials

### 3 layered cotton blend fabric

high density very stable hard-wearing especially skin friendly warming latex free

Size table  
Ellcare classic

kids

Adults

Size	XXS	XS	S	M	L
Circumf. (U) in cm	< 19	19 - 24	25 - 32	32 - 39	39 - 45
Length (L) in cm	35	38	45	51	60



# Ellcare flex

Item no. 322

## Elbow orthosis with incorporated aluminium splint, adjustable

Strong but bendable aluminium splint can be shaped as needed. Reverse drawing velcro straps for stabilisation of sides and adjusting the circumference.

Colours as shown in the picture. Children's sizes with velcro strap in different colours (non-eligible).



Children



Adults

### Indication

- After surgery
- Arthritis and arthrosis of elbow joint
- Instability of the elbow joint
- Posttraumatic and postoperative soft tissue injuries in and around the elbow
- Joint effusion and swelling

### Effectiveness

- Elbow stabilisation

### Measuring



**U = Circumference upper arm**  
Measure as shown

### Materials

**3 layered cotton blend fabric**  
high density  
very stable  
hard-wearing  
especially skin friendly  
warming  
latex free

Size table  
Ellcare flex



Adults

Size	XXS	XS	S	M	L
Circumf. (U) in cm	< 19	19 - 24	25 - 32	32 - 39	39 - 45
Length (L) in cm	35	38	39	42	48

# Epicare

Item no. 311

## Epicondylitis brace with incorporated pad for treatment of tennis elbow

This brace helps to relieve pain induced by muscle and tendon irritations or inflammations known as, for example, tennis elbow. The incorporated pad relieves the tendon on the spot.

Colour as shown in the picture.



### Indication

- Epicondylitis humeri (chronic or acute)
- Tennis elbow, golfers' elbow, mouse elbow

### Effectiveness

- pain-relieving
- tension-relieving
- strain-relieving

### Measuring



#### Measuring point

Circumference of the forearm, approx. 5–6 cm below the elbow

### Materials

#### Spacer fabric

#### 3 layered blend fabric

hygienic and stain resistant  
very hard-wearing  
very skin-friendly  
breathable  
sweat and heat regulating  
latex-free

### Size table Epicare

Adults

Size	S	M	L	XL
Circumf. (U) in cm	18 - 23	23 - 28	28 - 31	31 - 36



# HAND



Adults



## Carpfix plus

Hand and wrist support for children with volar aluminium splint

p. 28



## Manucare plus

Active bandage to support the wrist joint

p. 29



## Carpus classic plus

Hand and wrist support with volar aluminium splint

p. 30



## Carpus stabil plus

Hand and wrist support with volar aluminium splint for adults

p. 31



## Polli classic

Hand and thumb support with volar aluminium splint

p. 32



## Polli stabil plus

Hand and thumb support with volar aluminium splint

p. 33





# Carpfix plus

Item no. 331

## Hand and wrist support for children with volar aluminium splint

The hand and wrist support with integrated aluminium splint features velcro straps for additional fixation and reverse drawings which make it easier to put on.

Colours as shown in the picture. Children's sizes with velcro strap in different colours (non-eligible).



Children

### Indication

- Tendovaginitis
- Stabilisation after casts and surgery
- Chronic, posttraumatic and postoperative soft tissue conditions
- rheumatoid disorders
- Conditions after wrist joint distorsion
- Conditions after fractures (carpus incl. distal radius)
- Radial paresis
- Instability of the wrist joint
- Wrist arthrosis
- Carpal tunnel syndrome

### Effectiveness

- Stabilisation of the wrist joint with full mobility of the fingers at the same time

### Measuring



**U = Circumference wrist joint**

Measure as shown

### Materials

#### Spacer fabric

#### 3 layered blend fabric

hygienic and stain resistant  
very hard-wearing  
very skin-friendly  
breathable  
sweat and heat regulating

### Size table Carpfix plus



Size	XXXS	XXS	XS
Circumf. (U) in cm	< 11	11 - 13	13 - 15
Length (L) in cm	15	16	19
Side	left or right		



# Manucare plus

Item no. 332

**Active bandage to support the wrist joint with two incorporated pads, a volar synthetic splint as well as an elastic wrist joint strap**

The two incorporated pads support and relieve wrist joint as well as nerves and vessels. They also accelerate resorption of oedema and bruises. The volar synthetic splint stabilises the wrist joint and is removable for washing purposes.

The elastic strap can be positioned freely to the patients' needs and offers additional stability in the wrist joint area.

Excellent wear with highly elastic fabric, no slipping

Colour as shown in the picture.

**kids** + Adults




## Indication

- Tendovaginitis
- Arthrosis / arthritis
- posttraumatic
- Irritations (chronic, postoperative or posttraumatic)
- Conditions after wrist joint distorsion

## Effectiveness

- stabilising
- strain-relieving
- pain-relieving
- massaging
- supports healing

## Measuring



**U = Circumference wrist joint**

Measure as shown

## Materials

**3 layered synthetic blend fabric**

breathable

elastic

skin friendly

easy to clean

Pads: 100% silicon

Size table  
Manucare plus

**kids**

Adults

Size	XS	S	M	L	XL	XXL
Circumf. (U) in cm	14 - 15	15 - 16	16 - 17	17 - 18	18 - 19	19 - 21
Side	left or right	left or right				

# Carpus classic plus

Item no. 342

## Hand and wrist support with volar aluminium splint

The hand and wrist support with integrated aluminium splint features velcro straps for additional fixation and reverse drawings which make it easier to put on.

Colours as shown in the picture. Children's sizes with velcro strap in different colours (non eligible).

kids



Children



Adults


### Indication

- Tendovaginitis
- Stabilisation after casts and surgery
- Chronic, posttraumatic and postoperative soft tissue conditions
- rheumatoid disorders
- Conditions after wrist joint distorsion
- Conditions after fractures (carpus incl. distal radius)
- Radial paresis
- Instability of the wrist joint
- Wrist arthrosis
- Carpal tunnel syndrome

### Effectiveness

- Stabilisation of the wrist joint with full mobility of the fingers at the same time

### Measuring



**U = Circumference wrist joint**

Measure as shown

### Materials

#### Spacer fabric

#### 3 layered blend fabric

hygienic and stain resistant  
hard-wearing  
skin-friendly  
breathable  
sweat and heat regulating

### Size table Carpus classic plus

kids

Adults

Size	XXS	XS	S	M	L	XL
<b>Circumf. (U) in cm</b>	11 - 13	13 - 15	15 - 17	17 - 19	19 - 21	> 21
<b>Length (L) in cm</b>	19		19			
<b>Side</b>	left or right		left or right			

# Carpus stabil plus

Item no. 344

## Hand and wrist support for adults with volar aluminium splint

The hand and wrist support with integrated aluminium splint features velcro straps for additional fixation and reverse drawings which make it easier to put on.

Colour as shown in the picture.



Adults

### Indication

- Tendovaginitis
- Stabilisation after casts and surgery
- Chronic, posttraumatic and postoperative soft tissue conditions
- rheumatoid disorders
- Conditions after wrist joint distorsion
- Conditions after fractures (carpus incl. distal radius)
- Radial paresis
- Instability of the wrist joint
- Wrist arthrosis
- Carpal tunnel syndrome

### Effectiveness

- Stabilisation of the wrist joint with full mobility of the fingers at the same time

### Measuring



**U = Circumference wrist joint**

Measure as shown

### Materials

**Spacer fabric**  
**3 layered blend fabric**  
 hygienic and stain resistant  
 hard-wearing  
 skin-friendly  
 breathable  
 sweat and heat regulating

### Size table

#### Carpus stabil plus

Adults

Size	S	M	L	XL
<b>Circumf. (U) in cm</b>	15 - 17	17 - 19	19 - 21	> 21
<b>Length (L) in cm</b>	26 cm			
<b>Side</b>	left or right			

# Polli classic

Item no. 351

## Hand and thumb support with volar aluminium splint

Ideal fixation through two splints which can be adjusted according to anatomy. An additional velcro strap provides for perfect fitting.

Colours as shown in the picture. Children's sizes with velcro strap in different colours (non eligible).

kids



Children



Adults


### Indication

- Instability (skier's thumb)
- Arthrosis of the carpometacarpal joint
- Chronic, posttraumatic and postoperative irritations of the carpometacarpal joint

### Effectiveness

- Immobilisation of the metacarpophalangeal and carpometacarpal joint

### Measuring



**U = Circumference wrist joint**

Measure as shown

### Materials

**Spacer fabric**  
**3 layered blend fabric**  
 hygienic and stain resistant  
 hard-wearing  
 skin-friendly  
 breathable  
 sweat and heat regulating

Size table  
 Polli classic

kids

Adults

Size	XS	S	M	L	XL
Circumf. (U) in cm	13 - 15	15 - 17	17 - 19	19 - 21	> 21
Length (L) in cm	20	20			
Side	left or right	left or right			

# Polli stabil plus

Item no. 353

## Hand and thumb support with volar aluminium splint

Ideal fixation through two splints which can be adjusted according to anatomy. Additional velcro strap on the thumb.

Colour as shown in the picture.



Adults

### Indication

- Tendovaginitis
- Stabilisation after casts and surgery
- Chronic, posttraumatic and postoperative soft tissue conditions
- rheumatoid disorders
- Conditions after wrist joint distorsion
- Conditions after fractures (carpus incl. distal radius)
- Radial paresis
- Chronic, posttraumatic and postoperative irritations of the carpometacarpal joint
- Wrist arthrosis
- Carpal tunnel syndrome

### Effectiveness

- Stabilisation and immobilisation of the wrist joint with fixation of the metacarpophalangeal and carpometacarpal joint
- also usable as plaster cast replacement

### Measuring



**U = Circumference wrist joint**

Measure as shown

### Materials

**Spacer fabric**  
**3 layered blend fabric**  
 hygienic and stain resistant  
 very hard-wearing  
 Very skin-friendly  
 breathable  
 sweat and heat regulating

### Size table

#### Polli stabil plus

Adults










Size	S	M	L	XL
Circumf. (U) in cm	15 - 17	17 - 19	19 - 21	> 21
Length (L) in cm	26			
Side	left or right			





# FINGER

**kids** Adults

	<p><b>Pollifix</b> Thumb support with integrated aluminium splint p. 36</p>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Fingerfix</b> Strap for fixing two fingers or two toes to each other p. 37</p>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Stack finger splint</b> Finger splint for immobilising the DIP joints p. 38</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Stack finger splint - Starter Kit</b> All sizes in a box as basic equipment p. 38</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Stack finger splint for button hole deformity</b> Finger splint for immobilising the PIP joints p. 38</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Stack finger splint for button hole deformity - Starter Kit</b> All sizes in a box as basic equipment p. 38</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Fingersplint U shape</b> Aluminium splint with soft foaming p. 39</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Frog finger splint</b> Aluminium splint with soft foaming p. 40</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Baseball finger splint</b> Aluminium splint with soft foaming p. 41</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>


# Pollifix

Item no. 361

## Thumb support with integrated aluminium splint

Ideal fixation through two splints which can be adjusted according to anatomy. An additional velcro strap provides for perfect fitting.

Colours as shown in the picture. Children's sizes with velcro strap in different colours (non-eligible).



Children



Adults

### Indication

- Instability
- Arthrosis of the carpometacarpal joint
- Chronic, posttraumatic and postoperative irritations of the carpometacarpal joint

### Effectiveness

- Stabilisation and immobilisation of the metacarpophalangeal and carpometacarpal joint

### Measuring



**U = Circumference wrist joint**

Measure as shown

### Materials

**3 layered cotton blend fabric**  
 high density  
 very stable  
 hard-wearing  
 especially skin friendly  
 warming  
 latex free

### Size table Pollifix



Adults

Size	XXS	XS	S	M	L
Circumf. (U) in cm	< 11	11 - 13	13 - 15	15 - 17	17 - 23
Side	left or right			left or right	


# Fingerfix

Item No. 387


### Strap for fixing two fingers or two toes to each other

The loop design facilitates putting it on. Anti-slip foam holds the strap in place.

Colour as shown.



Fingerfix



Fingerfix narrow

### Indication

- Sprains
- Small fractures
- Jammed fingers

### Effectiveness

- Protection and stabilisation of finger or toe

### Measuring

**Fits all sizes**

### Materials

**Nylon, lycra, foam**

### Size table Fingerfix



Size	narrow	wide
Width in mm	13	26

## Stack finger splint

Item No. 385

## Stack finger splint for button hole deformity

Item No. 386

### Stack finger splint - Starter Kit

Item No. 385-5-10 or 385-5-20

Contains 8 pieces of all available sizes  
= 64 pieces



### Stack finger splint for button hole deformity - Starter Kit

Item No. 386-5-10 or 386-5-20

Contains 8 pieces of all available sizes = 56 pieces.



The finger splint made of perforated plastic material is available in two colours (clear and skin-coloured) and usable for all fingers (sizes 1 to 7) Both colours available.

#### Indication

- after fracture, luxation or ligament injury of the finger or thumb
- for protection from traumatic accidents
- extensor tendon ruptures
- Finger tip injuries (item no. 385 only)
- Nail bed injuries (item no. 385 only)

#### Effectiveness

- Immobilisation of DIP joints (item no. 385 only)
- Immobilisation of PIP joints (item no. 386 only)

#### Measuring

Item no. 385:

#### Circumference DIP joint

Item no. 386:

#### Circumference PIP joint

L = Product length in cm

#### Materials

**Unelastic plastic material**  
supporting shell material  
stabilising  
skin friendly  
easy to clean  
latex free

#### Size table Stack finger splint

Adults Item No. 385 | DIP joints

Adults Item No. 386 | PIP joint

Size	1	2	3	4	5	5.5	6	7	1	2	3	4	5	6	7
Length (L) in cm	5	5.3	5.5	6.2	6	6	6.5	6.2	5	5	5	5	5.4	5.4	5.7
Circumf. (U) in	4.5	5	5	6	6.5	7	7	7.8	4.5	5	5.5	6.5	7	7.5	7.9
Colour	clear or skin-coloured								clear or skin-coloured						



# Finger splint U shape / Finger splint universal

Item No. 381 / 381-2-0

**Aluminium splint with soft foaming**

This aluminium splint with soft foaming in U shape is ideal for quick and comfortable stabilisation of the finger after injury. It can be used for all fingers.

The finger splint U shape is also available as universal splint with a length of 50 cm. The splint can be cut and formed into shape as needed.

Colour as shown in the picture.



Adults

**Indication**

- Stabilisation of the interphalangeal joint after injuries of the front finger joint
- Fractures of the distal phalanx
- Fractures of the finger nail bed
- Severe contusions and sprains
- Injuries of the finger tip

**Effectiveness**

- Finger stabilisation

**Measuring**

L = finger length

**Materials**

**Aluminium with soft foaming latex-free**

Size table  
Finger splint U shape

Adults

Size	S	M	L	XL
Length (L) in cm	3.2	5.1	7.3	8.9

## Frog finger splint

Item no. 384

### Finger splint made of aluminium with soft foaming

The finger splint is ideal for quick and comfortable stabilisation of the finger after injuries. It can be used for all fingers.

Colour as shown in the picture.



Adults

### Indication

- Stabilisation of the interphalangeal joint after injuries of the front finger joint
- Fractures of the distal phalanx
- Fractures of the finger nail bed
- Severe contusions and sprains
- Injuries of the finger tip

### Effectiveness

- Stabilisation of the finger

### Measuring

L = finger length

### Materials

**Aluminium with soft foaming latex-free**

### Size table Frog finger splint

Adults

Size	S	M	L
Length (L) in cm	5.7	7.6	8.3

# Baseball finger splint

Item no. 383

## Finger splint made of aluminium with soft foaming

The finger splint is ideal for quick and comfortable stabilisation of the finger after injuries. It can be used for all fingers.

Colour as shown in the picture.



Adults

### Indication

- Stabilisation of the interpalangeal joint after injuries of the front finger joint
- Fractures of the distal phalanx
- Fractures of the finger nail bed
- Severe contusions and sprains
- Injuries of the finger tip

### Effectiveness

- Stabilisation of the finger

### Measuring

L = finger length

### Materials

Aluminium with soft foaming latex-free

### Size table

#### Baseball finger

Adults


Size	S	M	L
Length (L) in cm	8.9	10	11.4



# SPINE




Adults



### Lumbocare basic

Active support without pad for the muscular support of the lumbar vertebrae

p. 44



### Lumbocare basic plus

Active support with incorporated pad for the muscular support of the lumbar vertebrae

p. 45





# Lumbocare basic

Item no. 411

## Active support without pad for the muscular support of the lumbar vertebrae

This support for the lumbar spine acts supporting and stabilising.

Excellent wear with highly elastic fabric, no slipping

Colour as shown in the picture.



### Indication

- Static muscular insufficiency
- Lumbar spine syndrome
- Postoperative partial fixation

### Effectiveness

- stabilising
- strain-relieving
- pain-relieving
- massaging
- supports healing
- warming

### Measuring



**Measuring point**  
Circumference of the waist

### Materials

**Three layered synthetic blend fabric**  
skin friendly material  
latex-free

Size table  
Lumbocare basic



Adults

Size	XS	S	M	L	XL
Circumf. (U) in cm	70 - 80	80 - 90	90 - 100	100 - 110	110 - 120

# Lumbocare basic plus

Item no. 412

### Active support with incorporated pad for the muscular support of the lumbar vertebrae

This support for the lumbar spine acts supporting and stabilising. The triangular pad features massaging pimples to relieve tension and pain.

Excellent wear with highly elastic fabric, no slipping

Colour as shown in the picture.



#### Indication

- Static muscular insufficiency
- Lumbar spine syndrome
- Postoperative partial fixation

#### Effectiveness

- stabilising
- strain-relieving
- pain-relieving
- massaging
- supports healing
- warming

#### Measuring



**Measuring point**  
Circumference of the waist

#### Materials

**3 layered synthetic blend fabric**  
Pads: Silicon  
skin friendly material  
latex-free

Size table  
Lumbocare basic plus



Adults

Size	XS	S	M	L	XL
Circumf. (U) in cm	70 - 80	80 - 90	90 - 100	100 - 110	110 - 120



# KNEE

**kids** Adults

**kids** Adults



## Inocare knee splint

Knee orthosis with patella opening and lateral aluminium splints

p. 48



## Genocare plus

Active support for stabilisation and relief of the knee joint

p. 54



## Helpcare

Knee bandage with patella ring, patella opening and simple joint

p. 49



## Genocare ROM

Knee orthosis with aluminium splint and adjustable joint

p. 55



## Helpcare super

Patella bandage with patella ring, patella opening and cross fixation

p. 50



## Fixcare 0° / 10° / 20°

Three-part knee orthosis for immobilisation

p. 56



## Provicare

Knee bandage with patella ring, patella opening and two-axis joint

p. 51



## Patellacare

Patella strap for a targeted relief of the patella tendon

p. 57



## Provicare plus

Knee bandage with patella opening and two-axis joint for adults

p. 52



## Fourcare basic

Hard frame knee brace

p. 58



## Provicare senior

Knee orthosis with patella opening and two-axis joint for the elderly

p. 53



## Fourcare classic

Four point knee brace

p. 59



# Inocare knee splint

Item no. 525

## Knee orthosis for children with patella opening, lateral aluminium splints with physiological joint

In order to set the angle an Allen wrench is needed.  
You can set the following angles: 0°, 5°, 20°, 30° and 45° (extension) as well as 0°, 20°, 45°, 60° and 90° (flexion).

Colours as shown in the picture. Velcro strap in different colours per size (non-eligible).



### Indication

- Pre- and postoperative immobilisation
- Postoperative treatment in knee surgery
- Prevention of tibia subluxation
- After anterior cruciate injury or reconstruction
- For secure stabilisation injured or surgically treated ligaments

### Effectiveness

- Lateral and medial stabilisation of the knee
- Also usable as plaster cast replacement

### Measuring



**U = circumference knee centre**  
(measuring above patella)

### Materials

**Spacer fabric**  
**3 layered blend fabric**  
hygienic and stain resistant  
very hard-wearing  
very skin-friendly  
breathable  
sweat and heat regulating

### Size table

#### Inocare knee splint



Size	XXS	XS
Circumf. (U) in cm	27 - 30	30 - 33
Length (L) in cm	29	



# Helpcare

Item no. 541

## Knee bandage with patella ring, patella opening, lateral plastic splints and simple joint

This knee bandage stabilises the knee while allowing for full mobility of the joint. Four velcro fasteners openable towards the front facilitate putting on and off. Two flexible and adjustable all-around straps provide for additional stabilisation.

Colours as shown in the picture. Velcro straps in different colours per size (non-eligible).



Children

### Indication

- Oedema and swelling
- Insufficiency of knee joint muscles
- Patellofemoral pain syndrome
- Arthritis with ligament insufficiencies
- Gonarthrosis
- Functional protection after reconstruction of injured ligament structures

### Effectiveness

- Lateral and medial stabilisation of the knee
- Compression and proprioception effect

### Measuring



**U = circumference knee centre**  
(measuring above patella)

### Materials

**Spacer fabric**  
**3 layered blend fabric**  
 hygienic and stain resistant  
 very hard-wearing  
 very skin-friendly  
 breathable  
 sweat and heat regulating

Size table  
Helpcare



Size	XXXS	XXS	XS	S	M	L
Circumf. (U) in cm	21 - 23	23 - 25	25 - 27	27 - 29	29 - 31	31 - 33
Length (L) in cm	approx. 23					

# Helpcare super

Item no. 542

## Patella bandage with patella ring, patella opening and cross fixation

Patella bandage with patella ring for universal use. The bandage is closed all around. The cross fixation provides for secure fitting.

Colours as shown in the picture. Velcro straps in different colours per size (non-eligible).



### Indication

- Oedema and swelling
- Insufficiency of knee joint muscles
- Patellofemoral pain syndrome
- Arthritis with ligament insufficiencies
- Gonarthrosis
- Functional protection after reconstruction of injured ligament structures

### Effectiveness

- Compression and proprioception effect

### Measuring



**U = circumference knee centre**  
(measuring above patella)

### Materials

**Spacer fabric**  
**3 layered blend fabric**  
hygienic and stain resistant very hard-wearing  
very skin-friendly  
breathable  
sweat and heat regulating  
latex-free

### Size table

#### Helpcare super



Size	XXS	XS	S	M	L
Circumf. (U) in cm	23 - 25	25 - 27	27 - 29	29 - 31	31 - 33
Length (L) in cm	approx. 23				

# Provicare

Item no. 531

## Knee bandage with patella ring, patella opening, lateral plastic splints and two-axis joint

Two tight velcro fasteners for adjusting the circumference and on the back of the leg provide for ideal fitting. Two all-around straps create additional stability.

The Provicare can be used as rehabilitation splint or as sports splint.

Colours as shown in the picture. Velcro straps in different colours per size (non-eligible).



Children

### Indication

- pre- and postoperative immobilisation
- postoperative treatment in knee surgery
- prevention of tibia subluxation
- after anterior cruciate injury or reconstruction
- for secure stabilisation injured or surgically

### Effectiveness

- Lateral and medial stabilisation of the knee
- Compression and proprioception effect

### Measuring



**U = circumference knee centre**  
(measuring above patella)

### Materials

**Spacer fabric**  
**3 layered blend fabric**  
hygienic and stain resistant  
very hard-wearing  
very skin-friendly  
breathable  
sweat and heat regulating  
latex-free

Size table  
Provicare



Size	XXXS	XXS	XS	S	M	L
Circumf. (U) in cm	19 - 21	21 - 23	24 - 26	27 - 29	29 - 31	31 - 33
Length (L) in cm	23					

# Provicare plus


Item no. 533

## Knee bandage with patella ring, patella opening, lateral aluminium splints and two-axis joint

This bandage has three velcro fasteners above the patella to provide for ideal adjustment of the circumference. Four all-around straps offer additional fitting. In contrast to the Inocare knee splint the joint cannot be set to a certain angle. The all-around straps help to adjust the bandage ideally.

The Provicare plus can be used as rehabilitation splint or as sports splint.

Colours as shown in the picture. Children's sizes with velcro strap in different colours (non eligible).



Children



Adults

### Indication

- Oedema and swelling
- Insufficiency of knee joint muscles
- Postoperative, posttraumatic or rheumatic effusions
- Patellofemoral pain syndrome
- Arthritis with ligament insufficiencies
- Gonarthrosis
- Functional protection after reconstruction of injured ligament structures
- Medial rupture of collateral ligaments

### Effectiveness

- Lateral and medial stabilisation of the knee

### Measuring



**U = circumference knee centre**  
(measuring above patella)

### Materials

**Spacer fabric**  
**3 layered blend fabric**  
hygienic and stain resistant  
very hard-wearing  
very skin-friendly  
breathable  
sweat and heat regulating  
latex-free

### Size table Provicare plus

kids

Adults

Size	XXS	XS	S	M	L	XL	XXL	XXXL
Circumf. (U) in cm	30 - 32	32 - 34	33 - 35	35 - 37	37 - 39	39 - 42	42 - 45	45 - 49
Length (L) in cm	29		38					

# Provicare senior

Item no. 532

**Knee bandage with patella ring, patella opening, lateral plastic splints and two-axis joint, openable towards the front through zipper**

Two tight velcro fasteners for adjusting the circumference provide for ideal fitting. The Provicare senior can be used as rehabilitation splint or as sports splint.

The zipper facilitates putting the Provicare senior on and off, which is why it is especially suitable for elderly patients.

Colour as shown in the picture.



Adults

### Indication

- Oedema and swelling
- Insufficiency of knee joint muscles
- Pre- and postoperative immobilisation
- Patellofemoral pain syndrome
- Arthritis with ligament insufficiencies
- Gonarthrosis
- Functional protection after reconstruction of injured ligament structures
- Postoperative treatment in knee surgery
- Prevention of tibia subluxation
- After anterior cruciate injury or reconstruction
- For secure stabilisation injured or surgically treated ligaments

### Effectiveness

- Lateral and medial stabilisation of the knee

### Measuring



**U = circumference knee centre**  
(measuring above patella)

### Materials

**Spacer fabric**  
**3 layered blend fabric**  
 hygienic and stain resistant  
 very hard-wearing  
 very skin-friendly  
 breathable  
 sweat and heat regulating  
 latex-free

### Size table

#### Provicare senior

Adults

Size	S	M	L	XL	XXL	XXXL
<b>Circumf. (U) in cm</b>	33 - 35	35 - 37	37 - 39	39 - 42	42 - 45	45 - 48
<b>Length (L) in cm</b>	30	34				



# Genocare plus

Item no. 501

## Active support for stabilisation and relief of the knee joint with patella ring and two incorporated side bars

The active support helps to stabilise the joint and relief pain, thus supporting a faster mobilisation. With every move there is a massaging effect from the pad. Swellings and swelling tendencies are thus reduced. The patella is guided through the regular motion sequence by the pad.

Excellent wear with highly elastic fabric, no slipping.

Colour as shown in the picture.



### Indication

- Arthrosis / arthritis
- Fibromyalgia
- Ligamentosis
- Posttraumatic and postoperative irritations
- Relapsing joint effusion and swelling
- Patella instability
- Ligament laxity
- Gonarthrosis

### Effectiveness

- stabilising
- strain-relieving
- pain-relieving
- tension-relieving
- massaging
- supports healing

### Measuring



**U = Circumference knee**  
(measuring above patella)

### Materials

**3 layered synthetic blend fabric**  
breathable  
elastic  
skin friendly  
easy to clean  
Pads: 100% silicon  
latex-free

Size table  
Genocare plus



Adult

Size	XS	S	M	L	XL	XXL
Circumf. (U) in cm	25 - 28	28 - 31	31 - 34	34 - 37	37 - 40	40 - 44

# Genocare ROM


Item no. 521

## Knee orthosis with aluminium splint and adjustable joint

The setting of the extension or flexion limitation is done via pins in the joint scale (no further tools necessary). Four or six (depending on size) velcro fasteners for adjusting the circumference provide for ideal fitting.

The angles can be set in steps of 10 degrees without any further tools; 0° to 90° (Extension) and 0° to 120° (Flexion).

Colours as shown in the picture. Children's sizes with velcro strap in different colours (non eligible).



### Indication

- Pre- and postoperative immobilisation
- Postoperative treatment in knee surgery
- Prevention of tibia subluxation
- After anterior cruciate injury or reconstruction
- For secure stabilisation injured or surgically treated ligaments

### Effectiveness

- Lateral and medial stabilisation of the knee
- Also usable as plaster cast replacement

### Measuring



**U = circumference knee centre** (measuring above patella)

**L = overall length of orthosis**  
In relation to the overall length the joint is in the middle.

### Materials

**3 layered cotton blend fabric**  
high density very stable hard-wearing especially skin friendly warming latex free

### Size table Genocare ROM



Adults

Size	XXXS	XXS	XS	S	M	L	XL	XXL
Length (L) in cm	30	34	36	38	50	55	60	65
Circumf. (U) in cm	universal				universal			

## Fixcare 0°

Item no. 511

## Fixcare 10° *(no pediatric sizes)*

Item no. 514

## Fixcare 20°

Item no. 512

### Three-part knee orthosis for immobilisation with retral and lateral aluminium splint for children and adults

Knee orthosis for immobilisation. Due to the three-part design the orthosis is adjustable to the necessary circumference. The retral and lateral aluminium splints guarantee a high stability. Immobilisation, fixation and stabilisation of the knee. Knee pads available as accessories.

Colours as shown in the picture. Children's sizes with velcro strap in different colours (non-eligible).



Children



Adults

### Indication

- Pre- and postoperative immobilisation in knee surgery
- Patellofemoral pain syndrome
- Prevention of tibia subluxation
- Severe muscle injuries
- After anterior cruciate injury or reconstruction
- Condition after patella luxation
- Postoperative treatment after patellectomy
- Postoperative rehabilitation of injured or surgically treated ligaments (adults)
- Functional protection after reconstruction of injured ligament structures (adults)

### Effectiveness

- Immobilisation, fixation and stabilisation at 0°/10°/20°
- Also usable as plaster cast replacement

### Measuring

**L = L = overall length of orthosis**  
In relation to the overall length the joint is in the middle.

### Materials

**3 layered cotton blend fabric**  
high density  
very stable  
hard-wearing  
especially skin friendly  
warming  
latex free

Size table  
Fixcare 0° / 10° / 20°

**kids** (0° / 20° only)

Adults (Fixcare 10° available in M-XXL only)

Size	XXXS	XXS	XS	S	M	L	XL	XXL	XXXL
Length (L) in cm	26	30	40	40	50	55	60	65	70
Circ. (U) in cm	universal			universal					

# Patellacare

Item no. 502

## Patella strap for relief and stabilisation of the patella ligament

Patella bandage with incorporated pad. The bandage is put on in a way that the anatomically formed pad is situated below the patella in order to guide and support the patella ligament.

Colour as shown in the picture.



Adults

### Indication

- Pain in the patella area
- Tendinitis (ligament inflammation) of the patella ligament
- Patella tip syndrome
- Chondropathic patella
- Patellafemoral arthritis

### Effectiveness

- pain-relieving
- strain-relieving

### Measuring



#### Measuring point

Circumference below the patella

### Materials

**Spacer fabric**  
**3 layered blend fabric**  
 hygienic and stain resistant  
 very hard-wearing  
 very skin-friendly  
 breathable  
 sweat and heat regulating  
 latex-free

### Size table Patellacare

Adults

Size	S	M	L	XL
Circumf. (U) in cm	25 - 30	30 - 36	33 - 41	38 - 46

## Fourcare basic

Item no. 552-5-1xx

### Hard frame knee brace for adults

Knee brace made of robust aluminium; easy to adjust with pins (no wedges!) which are secured by a synthetic thread, NO tools necessary. The Fourcare basic is easy to handle as it is put on from one side only (frontside).

Colour as shown in the picture.

#### Possible settings:

Extension: 0°, 10°, 20, 40°, 60°, 90°

Flexion: 0°, 10°, 20, 40°, 60°, 90°, 150°



Adults

### Indication

- Complex knee instabilities, traumatic as well as degenerative
- Isolated or combined ligament instabilities of the knee joint
- After ligament reconstruction
- After meniscus reconstruction and resection
- Genu recurvatum
- Damages to or instabilities of the anterior and/or posterior crucial ligament (ACL, PCL, CL)
- Ligament ruptures or instabilities
- After patella fracture
- Long-term rehabilitation

### Effectiveness


- Prevention of tibial subluxation (anterior and posterior Drawer test)
- Reduction of inner rotation
- Limitation of extension and flexion
- postoperative and conservative ligament relief
- Immobilisation option

### Materials

#### 3 layered cotton blend fabric

high density  
very stable  
hard-wearing  
especially skin friendly  
warming  
latex free

### Measuring



**D = Diameter knee**



**U 1 = femoral circumf.**

approx. 10-15 cm above patella



**U 2 = lower leg circumf.**

approx. 10-15 cm below patella

**i Set of spare parts available**



### Size table

#### Fourcare basic

Adults

Size	XS	S	M	L	XL	XXL
<b>Diameter (D) in cm</b>	8.25 - 9.5	9.5 - 10.8	10.8 - 12.1	12.1 - 13.3	13.3 - 14.6	14.6 - 15.9
<b>Circumf. 1 (U) in cm</b>	35.6 - 42.4	42.4 - 50	50 - 54.6	54.6 - 61.6	61.6 - 68.6	68.6 - 75.6
<b>Circumf. 2 (U) in cm</b>	30 - 33.6	33.6 - 37.5	37.5 - 42.5	42.5 - 47.6	47.6 - 52.7	52.7 - 57.8
<b>Side</b>	left or right					



# Fourcare classic

Item no. 552-5-0xx

## 4 point hard frame knee brace for adults

Knee brace made of robust aluminium; easy to adjust with pins (no wedges!) which are secured by a synthetic thread, NO tools necessary. The Fourcare classic provides a high stability of the knee through opposite frame guidance.

Colour as shown in the picture.

### Possible settings:

Extension: 0°, 10°, 20, 40°, 60°, 90°

Flexion: 0°, 10°, 20, 40°, 60°, 90°, 150°



Adults

### Indication

- Complex knee instabilities, traumatic as well as degenerative
- Isolated or combined ligament instabilities of the knee joint
- After ligament reconstruction
- After meniscus reconstruction and resection
- Genu recurvatum
- Damages to or instabilities of the anterior and/or posterior crucial ligament (ACL, PCL, CL)
- Ligament ruptures or instabilities
- After patella fracture
- Long-term rehabilitation


### Effectiveness

- Prevention of tibial subluxation (anterior and posterior Drawer test)
- Reduction of inner rotation
- Limitation of extension and flexion
- postoperative and conservative ligament relief
- Immobilisation option

### Materials

**3 layered cotton blend fabric**  
 high density  
 very stable  
 hard-wearing  
 especially skin friendly  
 warming  
 latex free

### Measuring



**D = Diameter knee**



**U 1 = femoral Circumf.**  
 approx. 10-15 cm above patella



**U 2 = lower leg Circumf.**  
 approx. 10-15 cm below patella

**i Set of spare parts available**



### Size table

#### Fourcare classic

Adults

Size	XS	S	M	L	XL	XXL
<b>Diameter (D) in cm</b>	8.25 - 9.5	9.5 - 10.8	10.8 - 12.1	12.1 - 13.3	13.3 - 14.6	14.6 - 15.9
<b>Circumf. 1 (U) in cm</b>	35.6 - 42.4	42.4 - 50	50 - 54.6	54.6 - 61.6	61.6 - 68.6	68.6 - 75.6
<b>Circumf. 2 (U) in cm</b>	30 - 33.6	33.6 - 37.5	37.5 - 42.5	42.5 - 47.6	47.6 - 52.7	52.7 - 57.8
<b>Side</b>	left or right					













# LOWER LEG

**kids** Adults

	<p><b>Kids' Walker</b> Lower leg orthosis with lateral splints (without joint) for children p. 62</p>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
	<p><b>Walker classic / Walker classic low</b> Lower leg orthosis with lateral splints (without joint) for adults p. 63</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Walker One Piece</b> Lower leg orthosis with lateral splints without joint p. 64</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Walker ROM</b> Lower leg orthosis with lateral splints and joints p. 65</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Walker air</b> Lower leg orthosis with lateral splints without joint and air pad p. 66</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
	<p><b>Walker ROM air</b> Lower leg orthosis with lateral splints, joint and air pad p. 67</p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>

# Kids' Walker

Item no. 571

## Lower leg orthosis with lateral splints (without joint) for children

Synthetic inner shoe with cotton lining.

Easily adjustable to swellings through velcro fasteners, physiological gait pattern through rolloff sole.

Wounds can be inspected and treated easily, ease of use and comfortable to wear.

The inner part is available in red, green and blue.



Children

### Indication

- Arthropathia or neuropathia
- Trophic dysfunctions (ulcers, gangrene)
- Wound and fracture healing dysfunctions
- Postoperative treatment of ligament, soft tissue and tendon injuries
- Luxation fractures of the upper ankle joint
- Injuries of the forefoot or metatarsus as well as the tarsus
- Posttreatment of talus, calcaneus and medial malleolus fractures
- Stress fractures of tibia and fibula
- Stable distal tibia and fibula fractures
- For early functional therapies in Achilles tendon rupture
- After ligament refixation and peroneal tendon luxations
- Arthritis, e.g. infections arthritis, gout

### Size table Kids' Walker




Size	XXS	XS	S
Shoe size	20 - 23	24 - 26	27 - 30
Length (L) in cm	14	16.5	19
Height (H) in cm	17.8	20.3	22.9
Colour	please choose red, green, blue		

### Effectiveness

- stabilising

### Measuring



**L = Sole length**

Available in three colours: red, green and blue

### Materials

#### 3 layered cotton blend fabric

high density very stable hard-wearing especially skin friendly warming latex free

# Walker classic / Walker classic low

## Lower leg orthosis with lateral splints (without joint) for adults

Synthetic inner shoe with cotton lining.

Easily adjustable to swellings through velcro fasteners, physiological gait pattern through roll off sole.

Wounds can be inspected and treated easily, ease of use and comfortable to wear.

Colour as shown in the picture.



Adults

Walker classic



Adults

Walker classic low


### Indication

- Arthropathia or neuropathia
- Trophic dysfunctions (ulcers, gangrene)
- Wound and fracture healing dysfunctions
- Postoperative treatment of ligament, soft tissue and tendon injuries
- Luxation fractures of the upper ankle joint
- Injuries of the forefoot or metatarsus as well as the tarsus
- Posttreatment of talus, calcaneus and medial malleolus fractures
- Stress fractures of tibia and fibula
- Stable distal tibia and fibula fractures
- For early functional therapies in Achilles tendon rupture
- After ligament refixation and peroneal tendon luxations
- Arthritis, e.g. infections arthritis, gout

### Effectiveness

- stabilising

### Measuring



L = Sole length

### Materials

#### 3 layered cotton blend fabric

high density very  
stable hard-wearing  
especially skin friendly  
warming  
latex free

### Size table

Walker classic / Walker classic low Adults

Size	XS	S	M	L
Shoe size	30 - 33	34 - 38	38 - 42	43 - 46
Length (L) in cm	21.3	23.8	26.7	29.6
Height in cm, Walker classic	30	43	43	43
Height in cm, Walker classic low	24	24	25	25

# Walker One Piece

Item no. 572

## Lower leg orthosis with lateral splints without joint

Synthetic inner shoe with cotton lining.

Easily adjustable to swellings through velcro fasteners, physiological gait pattern through roll off sole.

Wounds can be inspected and treated easily, ease of use and comfortable to wear. Lateral splints and sole made of one piece for higher stability.

Colour as shown in the picture.



Adults


### Indication

- Arthropathia or neuropathia
- Trophic dysfunctions (ulcers, gangrene)
- Wound and fracture healing dysfunctions
- Postoperative treatment of ligament, soft tissue and tendon injuries
- Luxation fractures of the upper ankle joint
- Injuries of the forefoot or metatarsus as well as the tarsus
- Posttreatment of talus, calcaneus and medial malleolus fractures
- Stress fractures of tibia and fibula
- Stable distal tibia and fibula fractures
- For early functional therapies in Achilles tendon rupture
- After ligament refixation and peroneal tendon luxations
- Arthritis, e.g. infections arthritis, gout

### Effectiveness

- stabilising

### Measuring



L = Sole length

### Materials

#### 3 layered cotton blend fabric

high density very stable hard-wearing especially skin friendly warming latex free

### Size table

#### Walker One Piece

Adults

Size	S	M	L
Shoe size	34 - 38	38 - 42	43 - 46
Length (L) in cm	23.8	26.7	29.6
Height (H) in cm	43		



# Walker ROM

Item no. 575

## Lower leg orthosis with lateral splints and joints

Synthetic inner shoe with cotton lining.

Easily adjustable to swellings through velcro fasteners, physiological gait pattern through roll off sole.

Wounds can be inspected and treated easily, ease of use and comfortable to wear.

Colour as shown in the picture.



Adults


### Indication

- Arthropathia or neuropathia
- Trophic dysfunctions (ulcers, gangrene)
- Wound and fracture healing dysfunctions
- Postoperative treatment of ligament, soft tissue and tendon injuries
- Luxation fractures of the upper ankle joint
- Injuries of the forefoot or metatarsus as well as the tarsus
- Posttreatment of talus, calcaneus and medial malleolus fractures
- Stress fractures of tibia and fibula
- Stable distal tibia and fibula fractures
- For early functional therapies in Achilles tendon rupture
- After ligament refixation and peroneal tendon luxations
- Arthritis, e.g. infections arthritis, gout

### Effectiveness

- stabilising

### Measuring



L = Sole length

### Materials

**3 layered cotton blend fabric**  
 high density very stable hard-wearing especially skin friendly warming latex free

### Size table Walker ROM

Adults

Size	XS	S	M	L
Shoe size	30 - 33	34 - 38	38 - 42	43 - 46
Length (L) in cm	21.3	23.8	26.7	29.6
Height (H) in cm	30	43	43	43

## Walker air

Item no. 574

### Lower leg orthosis with lateral splints without joint and with air pad

Synthetic inner shoe with cotton lining. Additional pads available.

Easily adjustable to swellings through velcro fasteners, physiological gait pattern through rolloff sole.

Wounds can be inspected and treated easily, ease of use and comfortable to wear. Air padding provides perfect adjustability and compression to the lower leg.

Colour as shown in the picture.



Adults


### Indication

- Arthropathia or neuropathia
- Trophic dysfunctions (ulcers, gangrene)
- Wound and fracture healing dysfunctions
- Postoperative treatment of ligament, soft tissue and tendon injuries
- Luxation fractures of the upper ankle joint
- Injuries of the forefoot or metatarsus as well as the tarsus
- Posttreatment of talus, calcaneus and medial malleolus fractures
- Stress fractures of tibia and fibula
- Stable distal tibia and fibula fractures
- For early functional therapies in Achilles tendon rupture
- After ligament refixation and peroneal tendon luxations
- Arthritis, e.g. infections arthritis, gout

### Effectiveness

- stabilising

### Measuring



L = Sole length

### Materials

**3 layered cotton blend fabric**  
high density very  
stable hard-wearing  
especially skin friendly  
warming  
latex free

### Size table Walker air

Adult

Size	XS	S	M	L
Shoe size	30 - 33	34 - 38	38 - 42	43 - 46
Length (L) in cm	21.3	23.8	26.7	29.6
Height (H) in cm	30	43	43	43

# Walker ROM air

Item no. 576

## Lower leg orthosis with lateral splints, joint and air pad

Synthetic inner shoe with cotton lining. Easily adjustable to swellings through velcro fasteners, physiological gait pattern through roll off sole.

Wounds can be inspected and treated easily, ease of use and comfortable to wear. Air padding provides perfect adjustability and compression to the lower leg.

Colour as shown in the picture.



Adults


### Indication

- Arthropathia or neuropathia
- Trophic dysfunctions (ulcers, gangrene)
- Wound and fracture healing dysfunctions
- Postoperative treatment of ligament, soft tissue and tendon injuries
- Luxation fractures of the upper ankle joint
- Injuries of the forefoot or metatarsus as well as the tarsus
- Posttreatment of talus, calcaneus and medial malleolus fractures
- Stress fractures of tibia and fibula
- Stable distal tibia and fibula fractures
- For early functional therapies in Achilles tendon rupture
- After ligament refixation and peroneal tendon luxations
- Arthritis, e.g. infections arthritis, gout

### Effectiveness

- stabilising

### Measuring



L = Sole length

### Materials

**3 layered cotton blend fabric**  
 high density very stable hard-wearing especially skin friendly warming latex free

### Size table

#### Walker ROM air

Adults

Size	XS	S	M	L
Shoe size	30 - 33	34 - 38	38 - 42	43 - 46
Length (L) in cm	21.3	23.8	26.7	29.6
Height (H) in cm	30	43	43	43





# ANKLE JOINT

 Adults



## Castcare

Ankle joint orthosis, anatomically formed shell with shapeable foam material

p. 70



## Malleocare plus

Active support for the muscular ankle joint stabilisation featuring two silicon pads

p. 71



## Castcare

Item no. 561

### Ankle joint orthosis, anatomically formed shell with shapeable foam material

Shapeable foam material for ideal comfort of wear in the rehab phase. Flexible adjustment of circumference through velcro fasteners. No walking limitation. Also comfortable in shoes.

Colour as shown in the picture.



Children



Adults

### Indication

- Shorter plaster cast period
- Stabilisation of ankle joint
- Ligament injuries (fibular ligament rupture)
- Distorsions of the upper ankle joint
- Supination and/or lateral rotation trauma
- Stable lateral malleolus fractures

### Effectiveness

- Stabilisation and support of ankle joint area

### Measuring

#### Circumference universal

Two sizes: for children and adults

### Materials

#### 3 layered cotton blend fabric

high density very stable hard-wearing especially skin friendly warming latex free

Size table  
Castcare



Adults

Size

uni

uni



# Malleocare plus

Item no. 562

## Active support for the muscular ankle joint stabilisation featuring 2 silicon pads

The active support helps to stabilise the joint and relief pain, thus supporting a faster haematoma reduction. With every move the silicon pads provide a massaging effect, activating the blood circulation. Swelling predispositions or swellings are reduced and the healing accelerated.

Excellent wear with highly elastic fabric, no slipping

Colour as shown in the picture.



### Indication

- Irritations (postoperative and posttraumatic), e.g. After distortions
- Arthrosis / arthritis
- Relapsing joint effusion and swelling
- Ligament laxity
- Fibromyalgia
- After plaster cast treatments

### Effectiveness

- stabilising
- strain-relieving
- pain-relieving
- tension-relieving
- massaging
- supports healing

### Measuring



**Measuring point**  
Ankle circumference

### Materials

**3 layered synthetic blend fabric**  
breathable  
elastic skin  
friendly easy  
to clean  
Pads: 100% silicon  
latex-free

Size table  
Malleocare plus



Adults

Size	XS	S	M	L	XL	XXL
<b>Circumf. (U) in cm</b>	17 - 19	19 - 21	21 - 23	23 - 25	25 - 27	27 - 29
<b>Side</b>	left or right	left or right				



# FOOT

kids

Adults



### Universal kids' sandal / shoe

Children's support shoe with strong sole, open toes, sandal with open heel

p. 74



### Forefoot relief shoe

Forefoot relief shoe for relief and protection of the forefoot

p. 75



### Universal/Post surgery shoe

Universal shoe with strong sole

p. 76



### Plaster cast shoe

Shoe with strong sole and knuckle fixation

p. 77



## Universal kids' sandal

Item no. 582

## Universal kids' shoe (no picture)

Item no. 585

**Sandal / shoe for children with strong sole, slightly rounded and with semi-hard foam material, open toes. Closed heel (item no. 585) or open heel (item no. 582).**

Easy to use and comfortable to wear, even with larger dressings. Easy to adjust via velcro fasteners.

Both are available in red, green and blue.



Children


### Indication

- After Hallux valgus
- Fractures of the toes
- Injuries of the forefoot or metatarsus as well as the tarsus

### Effectiveness

- Forefoot relief after surgery and in wound care

### Measuring



**L = Length of sole**

### Materials

**3 layered cotton blend fabric**  
 high density very  
 stable hard-wearing  
 especially skin friendly  
 warming  
 latex free

### Size table

#### Universal kids' shoe



Size	XXS	XS	S	M
Shoe size	15 - 19	20 - 23	23 - 26	27 - 30
Length (L) in cm	11.4	14	16.5	19
Colour	please choose red, green, blue			

# Forefoot relief shoe

Item no. 581

## Shoe for relief and protection of the forefoot

One piece sole to protect the forefoot's underside. Pressure distribution to middle foot and heel. Can be widely opened to provide for a comfortable entrance. The width is adjustable through velcro strips.

Toe protection available as accessories.

Colour as shown in the picture.



Adults


### Indication

- After Hallux valgus
- Hammertoe or clawtoe surgery
- Healing of forefoot ulcerations

### Effectiveness

- Forefoot relief after surgery and in wound care

### Measuring



**L = Length of sole**

### Materials

- Nylon**
- elastic
- breathable (no sweating)
- skin friendly
- easy to clean
- latex free

### Size table Forefoot relief shoe

Adults

Size	XS	S	M	L	XL
Shoe size	33 - 36	36 - 39	39 - 42	42 - 45	45 - 48
Length (L) in cm	23	24.8	27.3	29	31

# Universal- / Post surgery shoe

Item no. 583

## Universal shoe with strong sole and semi-hard foam material

Easy to use and comfortable to wear, even with larger dressings. Easy to adjust via velcro fasteners.

Toe protection available as accessories.

Colour as shown in the picture.



Adults


### Indication

- After Hallux valgus
- Fractures of the toes
- Injuries of the forefoot or metatarsus as well as the tarsus

### Effectiveness

- Pressure relief, shock absorbing

### Measuring



**L = Length of sole**

### Materials

**Nylon**  
 elastic  
 breathable (no sweating)  
 skin friendly  
 easy to clean  
 latex free

### Size table

Universal / post surgery shoe Adults

Size	XS	S	M	L	XL
Shoe size	32 - 35	35 - 38	38 - 41	41 - 44	44 - 47
Length (L) in cm	22.5	24.4	26.3	28.6	30.5



## Plaster cast shoe

Item no. 584

### Shoe with strong sole, semi-hard foam material and knuckle fixation

Easy to use and comfortable to wear, even with larger dressings. A knuckle fixation provides for additional stabilisation. Easy to adjust via velcro fasteners.

Toe protection available as accessories.

Colour as shown in the picture.



Adults


### Indication

- After Hallux valgus
- Fractures of the toes
- Injuries of the forefoot or metatarsus as well as the tarsus
- Oedema

### Effectiveness

- Forefoot relief after surgery and in wound care

### Measuring



L = Length of sole

### Materials

**Nylon**  
 elastic  
 breathable (no sweating)  
 skin friendly  
 easy to clean  
 latex free

### Size table

#### Plaster cast shoe

Adults

Size	XS	S	M	L	XL
Shoe size	33 - 36	36 - 38	39 - 42	42 - 45	45 - 50
Length (L) in cm	22.8	24.2	26.8	28.5	33.4

# GENERAL TERMS AND CONDITIONS OF INOCARE MEDICAL GMBH

(scope of application: Germany, in German only)

## § 1 Allgemeines

1.1 Die nachfolgenden Bedingungen gelten für alle unsere Angebote, Verkäufe, Lieferungen und Leistungen und werden Inhalt des Vertrages. Sie gelten nicht, wenn unser Vertragspartner eine Privatperson ist und nicht beruflich oder gewerblich handelt. Sie gelten auch für alle künftigen Geschäftsbeziehungen, auch wenn sie nicht noch einmal ausdrücklich vereinbart werden.

1.2 Abweichenden oder ergänzenden Allgemeinen Geschäftsbedingungen des Käufers widersprechen wir hiermit ausdrücklich. Sie gelten auch dann nicht, wenn der Käufer sie seiner Bestellung oder sonstigen Erklärung zugrunde gelegt hat.

## § 2 Angebote und Aufträge

2.1 Unsere Angebote sowie Daten in unseren Preislisten sind frei bleibend, sofern sie nicht in schriftlicher Form als verbindlich bezeichnet sind. Ein wirksamer Vertrag kommt daher erst durch unsere Auftragsbestätigung oder die Auslieferung der Ware zustande.

2.2 Maßangaben, Gewichte, Abbildungen, Zeichnungen sowie andere Unterlagen, die zu unseren unverbindlichen Angeboten gehören bzw. dem Kunden zugänglich gemacht werden, bleiben in unserem Eigentum und sind nur annähernd maßgebend. Nur bei ausdrücklicher schriftlicher Bestätigung durch uns können sie verbindlicher Vertragsinhalt werden. Insbesondere begründen derartige Unterlagen keine Garantiehaftung.

## § 3 Zweifelhafte Zahlungsfähigkeit

3.1 Werden uns nach Vertragsschluss Umstände bekannt, die Zweifel an der Zahlungsfähigkeit des Käufers begründen, können wir weitere Lieferungen von einer Vorauszahlung der Ware durch den Käufer abhängig machen. Wir können dem Käufer für die Vorauszahlung der Ware eine angemessene Frist setzen und vom Vertrag zurücktreten, wenn die Vorauszahlung nicht fristgemäß bei uns eingeht; der Käufer kann statt der Vorauszahlung Sicherheit durch Bankbürgschaft leisten. Haben wir die Ware bereits geliefert, so wird der Kaufpreis ungeachtet vereinbarter Zahlungsfristen sofort ohne Abzug fällig.

3.2 Zweifel an der Zahlungsfähigkeit des Käufers sind unter anderem dann begründet, wenn ein Antrag auf Eröffnung des Insolvenzverfahrens über sein Vermögen gestellt wurde oder er Zahlungen an uns oder Dritte nicht pünktlich leistet.

## § 4 Preise

Unsere Preise verstehen sich netto in Euro ab Werk zuzüglich der jeweils gültigen gesetzlichen Umsatzsteuer.

## § 5 Lieferzeit

5.1 Alle genannten Liefertermine sind unverbindlich und gelten als nur annähernd vereinbart, soweit sie nicht von uns ausdrücklich als verbindlich bezeichnet worden sind. Bei unverbindlichen Lieferterminen gilt eine Lieferung innerhalb von drei Tagen nach der angegebenen Lieferzeit auf jeden Fall noch als rechtzeitig.

5.2 Wird uns die Leistung aufgrund höherer Gewalt oder aus anderen außergewöhnlichen und unverschuldeten Umständen ganz oder teilweise vorübergehend unmöglich oder erheblich erschwert, so verlängert sich die vereinbarte Lieferzeit um die Dauer des Leistungshindernisses. Gleiches gilt für eine gesetzliche oder vom Käufer gesetzte Frist für die Leistungserbringung, insbesondere für Nachfristen bei Verzug.

5.3 Bei einem etwaigen Lieferverzug, soweit er nicht auf Vorsatz oder grober Fahrlässigkeit beruht, sind Schadensersatzansprüche jeder Art ausgeschlossen.

## § 6 Versand / Versandkosten

6.1 Der Versand erfolgt auf Rechnung und Gefahr des Käufers, auch wenn frachtfreie Lieferung vereinbart ist. Wir sind nicht verpflichtet, für eine Transportversicherung zu sorgen.

6.2 Bei einem Nettoauftragswert von unter 250,00€ berechnen wir eine Versandkostenpauschale. Ab 250,00€ Nettoauftragswert je Bestellung liefern wir versandkostenfrei.

6.3 Bei speziellen Versandwünschen des Kunden (z.B. Nachnahme, Express, Samstagzustellung u. ä.) werden die uns entstehenden Kosten an den Kunden weiterberechnet.

6.4 Sofern nicht ausdrücklich schriftlich etwas anderes vereinbart wurde, sind wir zu Teillieferungen in zumutbarem Umfang berechtigt, die einzeln berechnet werden.

## § 7 Zahlung

7.1 Unsere Rechnungen werden generell elektronisch oder auf ausdrücklichen Wunsch des Kunden per Briefpost versendet und sind innerhalb von 20 Tagen ab Rechnungsdatum ohne Abzug zu zahlen. Bei Zahlung innerhalb von 10 Tagen nach Rechnungsdatum kann der Kunde 2% Skonto vom Rechnungsbetrag abziehen. Entscheidend ist hierbei der Zahlungseingang auf unserem Konto. Wir behalten uns vor, unberechtigten Skontoabzug nachzufordern.

7.2 Der Käufer kommt auch ohne eine Mahnung unsererseits in Verzug, wenn er den Kaufpreis nicht innerhalb von 3 Tagen nach Fälligkeit und Zugang der Rechnung oder einer gleichwertigen Zahlungsaufstellung zahlt.

7.3 Gerät der Käufer mit einer Zahlung in Verzug, werden seine sämtlichen Zahlungsverpflichtungen aus der Geschäftsverbindung mit uns sofort fällig. In diesem Fall sind wir berechtigt, von dem betreffenden Zeitpunkt an Zinsen in gesetzlich festgelegter Höhe zu verlangen. Der Nachweis eines höheren Schadens durch den Verkäufer bleibt vorbehalten.

7.4 Wir behalten uns vor, bei Zahlungsverzug Mahngebühren zu berechnen.

7.5 Der Käufer ist zur Aufrechnung, auch wenn Mängelrügen oder Gegenansprüche geltend gemacht werden, nur berechtigt, wenn die Gegenansprüche rechtskräftig festgestellt, von dem Verkäufer anerkannt wurden oder unstreitig sind. Zur Ausübung eines Zurückbehaltungsrechts ist der Käufer nur befugt, wenn sein Gegenanspruch auf demselben Kaufvertrag beruht.

## § 8 Retouren

8.1 Muster- und Auswahlendungen können unter Beifügung von Lieferschein und/oder Rechnung innerhalb der vereinbarten Frist zurück- gesendet werden. Sofern sich die Ware in einem einwandfreien und verkaufsfähigen Zustand befindet und die Rücksendung ausreichend frankiert wird, wird für derlei Retouren keine Kostenpauschale erhoben.

8.2 Für Retouren im Rahmen der Gewährleistung gilt §9 dieser AGB.

8.3 Retouren außerhalb der Gewährleistung werden grundsätzlich nur nach vorheriger Ankündigung und Abstimmung unter Beifügung von Lieferschein und/oder Rechnung entgegengenommen, wobei wir auch hier den einwandfreien Zustand voraussetzen. Ein generelles Verweigerungsrecht bei solchen Retouren behalten wir uns vor, insbesondere wenn der Zeitpunkt der Abnahme mehr als sechs Monate zurückliegt. Bei Gutschrift solcher Retouren erheben wir eine Wiedereinlagerungspauschale von mindestens 20% für Bearbeitung, Verpackung usw. Je nach Zustand kann dieser Satz höher ausfallen, bis zu 100%. Der Kunde hat das Recht, die Ware wieder zurückzuverlangen, falls für die Bearbeitung mehr als 20% einbehalten werden. Die hierfür entstehenden Versandkosten gehen zu Lasten des Kunden.

## § 9 Gewährleistung/Haftung

**9.1** Der Käufer hat die empfangene Ware auf Vollständigkeit, Transportschäden, offensichtliche Mängel, Beschaffenheit und deren Eigenschaften zu untersuchen. Gewährleistungsansprüche wegen offen- sichtlicher Mängel erlöschen, wenn sie nicht unverzüglich i. S. von §377 HGB uns gegenüber schriftlich geltend gemacht werden.

**9.2** Wir sind nicht zur Gewährleistung verpflichtet, wenn der Käufer einen offensichtlichen Mangel nicht rechtzeitig schriftlich gerügt hat. Soweit ein von uns zu vertretender Mangel an der Ware vorliegt und von dem Käufer rechtzeitig schriftlich gerügt wurde, sind wir – unter Ausschluss der Rechte des Käufers von dem Vertrag zurückzutreten oder den Kaufpreis herabzusetzen – zur Nacherfüllung verpflichtet, es sei denn, dass wir aufgrund der gesetzlichen Regelung zur Verweigerung der Nacherfüllung berechtigt sind. Der Käufer hat uns für jeden einzel- nen Mangel eine angemessene Frist zur Nacherfüllung zu gewähren.

**9.3** Die Nacherfüllung kann nach der Wahl des Käufers durch Beseitigung des Mangels oder Lieferung einer neuen Ware erfolgen. Wir sind berechtigt, die von dem Käufer gewählte Art der Nacherfüllung zu verweigern, wenn sie nur mit unverhältnismäßigen Kosten verbunden ist. Während der Nacherfüllung sind die Herabsetzung des Kaufpreises oder der Rücktritt vom Vertrag durch den Käufer ausgeschlossen. Eine Nachbesserung gilt mit dem zweiten vergeblichen Versuch als fehlgeschlagen. Ist die Nacherfüllung fehlgeschlagen oder hat der Verkäufer die Nacherfüllung insgesamt verweigert, kann der Käufer nach seiner Wahl Herabsetzung des Kaufpreises (Minderung) verlangen oder den Rücktritt vom Vertrag erklären.

**9.4** Schadensersatzansprüche zu den nachfolgenden Bedingungen wegen des Mangels kann der Käufer erst geltend machen, wenn die Nacherfüllung fehlgeschlagen ist oder die Nacherfüllung von uns verweigert wird. Das Recht des Käufers zur Geltendmachung von weitergehenden Schadensersatzansprüchen zu den nachfolgenden Bedingungen bleibt davon unberührt.

**9.5** Für vorsätzliche oder grob fahrlässige Pflichtverletzungen sowie für Schäden aus einer Verletzung des Lebens, des Körpers oder der Gesundheit haften wir uneingeschränkt nach den gesetzlichen Vorschriften. Im Übrigen haften wir nur, wenn die verletzte Vertragspflicht für das Erreichen des Vertragszwecks erkennbar von wesentlicher Bedeutung ist, und nur begrenzt bis zur Höhe des typischerweise vorhersehbaren Schadens.

**9.6** Die Haftungsbeschränkung nach Absatz 5 gilt entsprechend für andere als vertragliche Schadensersatzansprüche, insbesondere Ansprüche aus unerlaubter Handlung, mit Ausnahme der Ansprüche nach dem Produkthaftungsgesetz. Sie gilt ferner auch zugunsten unserer Angestellten, Arbeitnehmer, Mitarbeiter, Vertreter und Erfüllungsgehilfen.

**9.7** Soweit wir bezüglich der Ware oder Teile derselben eine Beschaffenheits- und/oder Haltbarkeitsgarantie abgegeben haben, haften wir auch im Rahmen dieser Garantie. Für Schäden, die auf dem Fehlen der garantierten Beschaffenheit oder Haltbarkeit beruhen, aber nicht unmittelbar an der Ware eintreten, haften wir allerdings nur dann, wenn das Risiko eines solchen Schadens ersichtlich von der Beschaffenheits- und Haltbarkeitsgarantie erfasst ist.

**9.8** Wir haften auch für Schäden, die durch einfache Fahrlässigkeit verursacht werden, soweit diese Fahrlässigkeit die Verletzung solcher Vertragspflichten betrifft, deren Einhaltung für die Erreichung des Vertragszwecks von besonderer Bedeutung ist (Kardinalpflichten). Wir haften jedoch nur, soweit die Schäden in typischer Weise mit dem Vertrag verbunden und vorhersehbar sind. Bei einfachen fahrlässigen Verletzungen nicht vertragswesentlicher Nebenpflichten haften wir im Übrigen nicht. Die in §7 enthaltenen Haftungsbeschränkungen gelten auch, soweit die Haftung für die gesetzlichen Vertreter, leitenden Angestellten und sonstigen Erfüllungsgehilfen des Verkäufers betroffen ist.

**9.9** Eine weitergehende Haftung ist ohne Rücksicht auf die Rechtsnatur des geltend gemachten Anspruchs ausgeschlossen. Soweit die Haftung des Verkäufers ausgeschlossen oder beschränkt ist, gilt dies auch für die persönliche Haftung seiner Angestellten, Arbeitnehmer, Mitarbeiter, Vertreter und Erfüllungsgehilfen.

## § 10 Eigentumsvorbehalt

**10.1** Wir behalten uns das Eigentum an der Ware vor (Vorbehaltsware), bis zum Eingang aller Zahlungen aus dem Kaufvertrag. Die gelieferten Waren gehen erst dann in das Eigentum des Käufers über, wenn dieser seine gesamten Verbindlichkeiten aus der Geschäftsverbindung einschließlich Nebenforderungen, Schadensersatzansprüchen und Einlösungen von Schecks und Wechseln erfüllt hat.

**10.2** Der Käufer hat uns von allen Zugriffen Dritter, insbesondere von Zwangsvollstreckungsmaßnahmen sowie sonstigen Beeinträchtigungen seines Eigentums unverzüglich schriftlich zu unterrichten. Der Käufer hat uns alle Schäden und Kosten zu ersetzen, die durch einen Verstoß gegen diese Verpflichtung und durch erforderliche Maßnahmen zum Schutz gegen Zugriffe Dritter entstehen.

**10.3** Kommt der Käufer seiner Zahlungsverpflichtung trotz einer Mahnung unsererseits nicht nach, so können wir die Herausgabe der noch in seinem Eigentum stehenden Vorbehaltsware ohne vorherige Fristsetzung verlangen. Die dabei anfallenden Transportkosten trägt der Käufer. In der Pfändung der Vorbehaltsware durch uns liegt stets ein Rücktritt vom Vertrag. Wir sind nach Rückbehalt der Vorbehaltsware zu deren Verwertung befugt. Der Verwertungserlös wird mit unseren offenen Forderungen aufgerechnet.

## § 11 Erfüllungsort

Der Erfüllungsort für Zahlungen ist der Sitz unseres Unternehmens, für unsere Warenlieferungen der Versandort.

## § 12 Datenverarbeitung

Der Käufer ist damit einverstanden, dass wir die im Zusammenhang mit der Geschäftsbeziehung erhaltenen Daten über den Käufer unter Beachtung des Bundesdatenschutzgesetzes für die Erfüllung eigener Geschäftszwecke verarbeiten, insbesondere speichern oder an eine Kreditutzschutzorganisation übermitteln, soweit dies im Rahmen der Zweckbestimmung des Vertrages erfolgt oder zur Wahrung unserer berechtigten Interessen erforderlich ist und kein Grund zu der Annahme besteht, dass das schutzwürdige Interesse des Käufers an dem Ausschluss der Verarbeitung, insbesondere der Übermittlung, dieser Daten überwiegt.

## § 13 Gerichtsstand und anzuwendendes Recht

Für das Vertragsverhältnis zwischen dem Käufer und uns gilt ausschließlich das Recht der Bundesrepublik Deutschland, auch wenn der Käufer seinen Wohn- oder Geschäftssitz im Ausland hat. Die Anwendung des einheitlichen Gesetzes über den internationalen Kauf beweglicher Sachen sowie des Gesetzes über den Abschluss von internationalen Kaufverträgen über bewegliche Sachen ist ausgeschlossen.



**inocare Medical GmbH**

Industriestraße 42  
74927 Eschelbronn, Germany

Tel.: +49 6226 7873 222  
Fax: +49 6226 7873 224

Mail:  
info@inocare.de  
www.inocare.de

